

Vials

Sampling and Storage Containers

Large Amber Screw Top Septum Bottles

Protect your sample from sunlight. Only glass and Teflon contact the contents.



DESCRIPTION	QTY.	CAT. NO.	PRICE
30mL (32 x 80mm)			
Bottles w/hole cap, Teflon/silicone septa	12	23228	
Bottles	12	23232	
Hole caps, no septa	100	27176	
Solid caps, polyseal liner	144	23238	
Septa, Teflon/silicone	144	23242-U	
60mL (38 x 95mm)			
Bottles w/hole cap, Teflon/silicone septa	12	23229	
Bottles	12	23233	
Hole caps, no septa	100	27176	
Solid caps, polyseal liner	144	23238	
Septa, Teflon/silicone	144	23242-U	
120mL (49 x 114mm)			
Bottles w/hole cap, Teflon/silicone septa	12	23230-U	
Bottles	12	23234	
Hole caps, no septa	144	23236	
Solid caps, polyseal liner	144	23240-U	
Septa, Teflon/silicone			
2 mil Teflon, 0.120"	100	27236	
10 mil Teflon, 0.100"	100	27237-U	
240mL (61 x 143mm)			
Bottles w/hole cap, Teflon/silicone septa	12	23231	
Bottles	12	23235	
Hole caps, no septa	144	23237	
Solid caps, polyseal liner	144	23241	
Septa, Teflon/silicone	144	23245-U	

Screw Top Mini Flasks

Heat resistant borosilicate glass flasks provide an inert seal when used with a hole cap and Teflon silicone septum (included) or a Mininert valve. Excellent for liquid-liquid extraction of trihalomethanes from drinking water. Also useful for sample collection, solution preparations, derivatization reactions, extractions, and many other applications.



DESCRIPTION	QTY.	CAT. NO.	PRICE
10mL	12	64715	
25mL	12	64716-U	
Screw caps w/hole, 20mm	12	27176	
Septa, PTFE/silicone	100	27177	
Mininert valves, 20mm	12	33303	



Crimpers and Decappers

DESCRIPTION	CAT. NO.	PRICE
8mm		
Hand Crimper	33272-U	
Pliers-type Decapper	33284	
Hand Decapper	33180-U	
11mm		
Hand Crimper	33195	
Hand Crimper, Adjustable	22313-U	
Pliers-type Decapper	33281	
Hand Decapper	33181-U	
13mm		
Hand Crimper	33279-U	
Pliers-type Decapper	33282	
Hand Decapper	33182-U	
20mm		
Hand Crimper	33280-U	
Hand Crimper, Adjustable	22316-U	
Pliers-type Decapper	33283	
Hand Decapper	33183	